

Technical Sheet

Carisolv®

RLS Global AB

Product Indication

Product type	Caries Removal Gel
Intended use/ purpose	Softens demineralized carious dentine prior to mechanical removal during chemo-mechanical minimal invasive treatment of patients diagnosed with caries.

Product properties

Origin	Sweden
Storage temp.	Keep the product refrigerated (2-8°C) in its original packaging during storage.
Product Specifications	Carisolv consists of two solutions provided in a double syringe with a closable cap. The solutions are mixed to form a gel at the time of use. Each package consists of 2 syringes packaged in aluminium pouches, 2 plungers, 10 static mixers and 10 intra-oral tips. The mixing of the solutions forms a clear gel, with a slight chlorine odour.
Ingredients	A clear, viscous solution consisting of water, carboxymethyl cellulose, sodium hypochlorite, sodium chloride, and amino acids (glutamic acid, leucine, lysine) with alkaline pH.
Sterilization / Desinfection	N/A
Shelf life	15 months
Regulatory status	EU: Medical Device Class I
Entry on the market	1998

Packaging

Packaging	2 pieces of Carisolv units are packaged together with 2 plungers, 10 static mixers, 10 intra-oral tips and an IFU in a shelf box.
-----------	---

Available Product line

Product	Type/Size/Amount	Art.-No.
---------	------------------	----------

Carisolv®	2 pcs twin syringes containing 2.5 ml gel, 10 pcs mixing elements and 10 pcs applicator tips For up to 8 treatments depending on size of the cavity.	10491
-----------	---	-------

Customer Contact

RLS Global AB
Customer contact
Tel.
+46 31-780 68 20
E-mail
customer@rlsglobal.se