

Retreatment & Repair



EVS540 - EVS545

S5 Retreatment Rotary Files

The removal of old filling materials is achieved quickly and easily. Simply apply the Crown Down technique using only 4 files (sizes 1 to 4).

- Markings of the S5 Retreatment files correspond to the programs of the S5 Endo motor.
- Files available in 23 mm length.

Sendoline
Perfect Endo



EJM005 - EJM008

Sterile ProTaper Retreatment

ProTaper retreatment files are designed to remove filling material before canal reshaping. The system consists of 3 instruments specifically designed for easy un-filling:

- D1 (16mm) – for coronal filling removal,
- D2 (18mm) – for middle filling removal
- D3 (22mm) – for apical filling removal.

Assorted pack also available. Pack of 6.

- Dentsply Sirona



EVD075

Biodentine

All-in-one, biocompatible and bioactive material for use wherever dentine is damaged.

- Based on unique Active Biosilicate Technology designed to treat damaged dentine both for restorative and endodontic indications.
- Fully biocompatible mineral material that reduces any risk of tissue response and is also a bioactive material promoting reactionary dentine genesis for pulp vitality preservation.
- For crown and root indications.

Box of 15 Capsules.

- Septodont



EVD071

Endosolv

Solvent for root canal sealers. Indications:

- Solvent for endodontic cements.
- This solvent is able to soften conventional zinc oxide eugenol cements and phenolic resin type sealers in the context of endodontic retreatment.
- Helps to remove zinc oxide eugenol-based and phenolic resin type root canal sealers.
- One single product instead of two.
- High softening potency.
- Reduces risks associated with use of power driven instruments during mechanical desobturation.



EVM012 / EVM015

Pro Root MTA Mineral Trioxide Aggregate

Root canal repair material.

- Water-based chemistry allows normal setting in the presence of moisture.
- Resistance to marginal leakage and reduction of bacterial migration.

Available in 2 dose and 5 dose packs.

- Dentsply Sirona

Call 0800 585 586 | www.ddgroup.com

TheEndoCollection

EVA000

MTA Capsules

Based on Mineral Trioxide Aggregate, MTA Caps is a proven endodontic repair material with great clinical results.

- Repairs root perforations and is used for apical filling in endodontic surgical procedures.
- Well known for effective pulpal protection following pulpotomies, particularly in children, MTA Caps is perfectly indicated as a treatment choice for apexification and apexogenesis procedures.

Contains 2 x 0.3g capsules.

- Acteon

