## **Surgical Procedure**

Luxators & Elevators



## **Luxator Periotome Titanium**

The Luxator Periotome Titanium offers enhanced durability. Extraction is performed with minimum tissue damage. Titanium coating - glides easier into the socket. The durable titanium coated tip does not require sharpening. The DE3 Titanium Dual Edge enables safer stepped insertion, reducing the risk of slipping. Available in 4 blade sizes: L3S 3mm Titanium Straight Blade, Grey L2S 2mm Titanium Straight Blade, Green DE3 3mm/1.5mm, Titanium Dual-Edge, Grey L3A 3mm, Titanium Angled Blade. Lilac





Contains 1 of each: L3C, L3S, L5C, L5S and a sharpening tool.

**DEE200** 

Luxator Periotome Kit L-K7

Includes: Luxator L2S, L3C, L3S, L3CA, L3IC, L5C, L5S and Luxator sharpening tool L10.

**DEE210** 

# LX Mechanical Periotome Starter Kit

Setting a new standard for trauma free extraction, the Luxator LX mechanical periotome offers:

- · Improved safety, precision and access
- Quicker extractions with minimal tissue damage, helps access hard to reach areas.
- The self directing tip allows the periotome tip to follow the root surface, reducing the risk of bone damage.
- A reciprocating tip enables a secure cutting of the fibers in a patient-friendly way.
- The Luxator LX mechanical periotome will only cut when under pressure, which means no risk of slipping.
- The titanium coated tip is a durable tip and remains sharp.

Contains: 1 Handpiece, 2 Periotome Tips Short 2mm & 3mm, 2 Periotome Tips Long 2mm & 3mm, tip holder & plunger

**DEE220** 







DEE000

#### **Premier Elevators**

Couplands Couplands No. 1 No. 3

DEE411 DEE413

Couplands Howarth No. 2 DEE350

A wider range of UnoDent Elevation Instruments are available online at ddgroup.com

#### **Elevators**

Warwick James Straight. Used to loosen the tooth from the periodontal ligament and ease extraction. **DEE511** 

### **Luxating Kit**

Contains 1 of each: EL3C, EL3S, EL5C, EL5S, Sharpening Stone.

DEE520

## **Hybrid Luxating Kit**

Contains 1 of each: EPTLGD, EPTLGM, EPTSMC, EPTSMD, EPTSMM, EPTSMS, EPTSSP, Sharpening Stone and IME12DIN86 cassette.

**DEE521** 

