

Obturation

Sendoline
Perfect Endo



S5

EVS580 / EVS582

Paper Points

Rolled to exact dimensions optimising the root canal treatment with its excellent absorbance ability. Available in Size 2 (30 Taper .06) & Size 3 (30 Taper .04). Pack 60.

EVS572-EVS573

Gutta Percha Points

Specifically designed to match the canal shape created with the S5 Rotary Files. Finish the preparation with either file No.2 or No.3 and simply select the S5 GP that corresponds to that file. It is precise, rapid and easy to use with only one Gutta Percha point required. Available in Size 2 (30 Taper .06) & Size 3 (30 Taper .04). Pack 60.

- Sendoline

EPK080-EPK085

Paper Points

Available in Small 2 to 3, Medium/Large 1 to 3 and Assorted. Pack 100

EQK330-EQK336

Gutta Percha Points

Available in Small 2 to 3, Medium/Large 1 to 3 and Assorted. Pack 50

TF Adaptive

TF Adaptive Gutta Percha and Paper Points are made specifically for matching TF Adaptive Files sizes and tapers with the same colour-coding.



Kerr
ENDODONTICS

EVE030 / EVE032

AH Plus Jet

AH Plus Jet sealer delivers added efficiency and precision; its preloaded sealer remains true to its proven, patented formula.

- With excellent biocompatibility, it is a viable solution to post-operative problems such as sensitivity and peri-apical inflammation.
- Low solubility makes it resistant to breaking down over time.
- Minimal shrinkage reduces gapping between the sealer and the canal wall.
- Should the need for removal arise, the sealer's thermoplastic properties make the task easier, whether using chemicals or a drill.

- Dentsply Sirona



WaveOne Gold

EPM470-EPM476

Sterile Paper Points

- Dedicated WaveOne Gold paper points.
- Colour-coded handles

Available in Small, Primary, Medium & Large. Pack of 180

EQM080-EQM084


Conform Fit Gutta Percha Points

Perfectly match the size of the canals prepared with WaveOne Gold files, letting you achieve rapid and precise warm 3D filling. Colour-coded handles provide a simple transition from instrumentation to obturation. Heated gutta-percha combined with the central carrier ensures excellent filling of the apex and lateral canals. Minimal risk of root fracture. Single insertion simplicity is more effective than lateral condensation. Less technique-sensitive than warm vertical compaction procedures.

Gutta-Percha Points designed for WaveOne Gold Files and for all warm or cold obturations

ISO Colour-Coded, Matches WaveOne Gold Files system, 29mm Length, Pack 60





EVK450
SendoSeal MTA

SendoSeal MTA is a root canal sealer designed to permanently seal the root canal. The sealer is bio-compatible and MTA-based, which means it sets by absorbing moisture in the root canal. With its high pH value (pH 10–12), SendoSeal MTA has an excellent antibacterial effect, which helps with the healing process.

- Simple Application – SendoSeal MTA root canal sealer comes premixed in a syringe to apply the material directly to the root canal with the thin cannula, which ensures smooth and simple handling. It has an outstanding flowability, making it easy to fill the entire root canal.
- Dimensional Stability and Hermetic Sealing – SendoSeal MTA root canal sealer is dimensionally stable and does not dissolve over time. The root canal remains sealed, protecting against the risk of any new infection.
- High Radiopacity – SendoSeal MTA has a high radiopacity of over 8 mm Al which makes it easy to distinguish the root filling and assess the result of the X-ray.

Box of 1x 3g syringe Sendoseal MTA
 + 20 x 24 gauge cannulas.

Sendoline
 Perfect Endo



EVK380 - EVK395

TubliSeal

TubliSeal EWT / TubliSeal Xpress / EWT Xpress

TubliSeal is a paste/paste zinc oxide eugenol root canal sealer. Also available TubliSeal Xpress in Automix syringes to help simplify setup, delivery and clean-up. This innovative automixing sealer adds efficiency, convenience and simplicity to a line of canal sealants that is already steadfast and dependable.

- TubliSeal is light in colour, is non-darkening and is radio-opaque.
- TubliSeal EWT offers a longer processing time.

- Kerr Endodontics

Kerr
 ENDODONTICS

EVK340 / EVK345

Sealapex

Sealapex Xpress

Sealapex is the original non-eugenol, calcium hydroxide polymeric root canal sealant.

- This formulation produces rapid healing and hard tissue formation.
- Sealapex Xpress formulation is available in Automix Syringes for better compliance on cross-infection prevention.

- Kerr Endodontics

Call 0800 585 586 | www.ddgroup.com

TheEndoCollection

EQS314

Elements Free

Cordless Continuous Wave Obturation

Designed for use with the warm vertical condensation technique, the Elements Free obturation system offers both downpack and backfill capabilities in a cordless design – giving dentists and endodontic specialists the freedom of movement to perform endodontic procedures anywhere without restrictions.

- Features a 360° one-touch activation ring to instantly heat the plugger and quickly cool upon deactivation.
- Digital temperature control for procedural accuracy.
- Downpack heats to 200°C in less than 0.5 seconds and with an adjustable heating range of 100 – 400°C.
- Swiss motor and gear drive for exceptional performance.
- Motorised extrusion for less hand fatigue
- Single-use cartridges with built-in finger grips for ease of use.
- Aerogel insulation surrounds the heating element for your comfort and your patient's safety.
- Perfect for Continuous Wave and single motion downtick obturation.





EVU005

UnoCal NS

Pre-mixed Calcium Hydroxide, Ca(OH)₂, dressing paste with Barium Sulphate for temporary canal filling. Features: Easy to remove with good water solubility, Excellent antibacterial effect, Superior radio-opacity, Strong bactericidal effect, High flowability, Excellent biocompatibility, Easy to handle (syringe paste type) Non-Setting Ca(OH)₂ Paste

Contains: 2g syringe, disposable tip, Endo-stop & Tip cap



EVU000

UnoSeal R

Resin-based sealer for the permanent obturation of root canal. Compatible with any heated or cold obturation technique. Features: Excellent biocompatibility, Superior radio-opacity, Tight sealing, Easy to mix paste-paste (1:1), Insoluble to tissue fluids, Excellent antibacterial effect, Lubricant during insertion of GP point

Contains: 1 x 10g Dual syringe (Base 5g, Catalyst 5g), Mixing tip, Mixing pad & Eco tip




LII005

Unocal LC

Light-curing Radiopaque Calcium Hydroxide Unocal/Unocal LC are self and light-curing radiopaque calcium hydroxide. Excellent adhesion – especially to composites. Excellent biocompatibility. High compressive strength and very stable. High resistance to amalgam condensation. Release of calcium hydroxide promotes secondary dentine creation.

Contains: 2 x 4g syringes and 20 x intra-oral tips



EVE036

AH Plus Bioceramic Sealer Starter Kit

AH Plus® Bioceramic Sealer is a calcium silicate-based root canal sealer. Tricalcium Silicate interacts with body fluids to release calcium and hydroxide ions that promote hydroxyapatite (HA) formation and supports the healing response.

It is available in a convenient, easy-to-use syringe with a small, 24-gauge cannula that can be placed directly into the canal.

AH Plus® Bioceramic Sealer is compatible with, and supplements, the most innovative Endodontic treatment solutions of Dentsply Sirona, e.g. WaveOne® Gold, ProTaper Ultimate™ and TruNatomy®.

Key Features & Benefits:

- A faster setting seale
- Dimensionally stable: low solubility of 0.11% with a high washout resistance to seal the root canal
- Safe and biocompatible: does not stimulate periodontal tissue.
- Better radiopacity
- Enhanced aesthetics: free of bismuth oxide.

Starter kit includes: 1 x 3g Syringe, 20 x Cannulas, 5 x Storage Caps and User Card



LII000

UnoSeal SC

Unocal Self-curing Radiopaque Calcium Hydroxide Unocal/Unocal LC are self and light-curing radiopaque calcium hydroxide. Excellent adhesion – especially to composites. Excellent biocompatibility. High compressive strength and very stable. High resistance to amalgam condensation. Release of calcium hydroxide promotes secondary dentine creation. Latex-free.

Contains: 1 x 13g Tube Base, 1 x 11g Tube Catalyst and Mixing Pad.



EVD080

BioRoot RCS



Permanent Root Canal Sealer.

- Outstanding adhesion to dentine and guttapercha points.
- Hydrophilic, continues the sealing process in the presence of moisture.
- Pure mineral formulation, will not stain teeth.
- Resin-free, made from pure calcium silicate and is monomer-free ensuring zero shrinkage.



- Great flowability, seals auxiliary canals.
- Uses cold single cone or cold lateral condensation.
- High 5mm Al radio-opacity for clear images on radiographs

