## **Surgical Procedure**

Instruments





## **Surgical Instruments**

Every Anaqa product is subjected to a rigorous four-step quality control process. This robust process is the reason why Anaqa are able to confidently offer industry-leading warranties on all their instruments.

- Slim beaks for better access and hold
- Stunning finish with high corrosion resistance
- Four stage QC process
- Optional titanium coating for premium finish and maximum corrosion resistance
- High quality medical grade steel
- Smooth joints for stress-free working
- Ergonomic, lightweight design
- Industry-leading warranties

#### **Luxation Instrument 5ST**

Black Titanium Nano Coating. DEE803

#### **Periotome**

Black Titanium Nano Coating Available in sizes 1, 2 & 3 DEE814-816

### Millers Bone File No. 52

Black Titanium Nano Coating Straight cut bone file used for final smoothing of bone.

### **Periotome Kit**

Used to sever the gingival fibres and periodontal ligament. Once a tooth is removed the Periotome can be used to mobilise the soft tissue. Minimises treatment of the alveolar bone and comes with 3 blades to facilitate access throughout the mouth. Straight blades for anterior teeth and premolars. Angled blade for posterior teeth.

Contains - One standard stainless steel handle for detachable tips, one periotome blade narrow, one periotome blade wide and one periotome blade angled all in standard stainless steel.

#### DEE818

Blade P1 Wide DEE820 Blade P2 Narrow DEE819 Blade P3 Angled DEE821 Handle DEE822





# Micro Surgical Basic Kit

The Kit includes:

IMS Infinity Cassette 1 DIN for 18 instruments, Black, High Definition Mirror, Mirror Handle, UNC probe 1-15mm, Micro Surgical Scalpel Handle SinusLine, Micro Tissue Plier, Periosteal Micro Surgery, Periosteal Trombelli, Micro Surgical Scissors Cortellini SuperCut. Needle holder

Castroviejo
Perma Sharp,
Suture plier,
Urban Cheek
Retractor,
Labanca
Malleable
Retractors,
Set of 3.

**DDH008** 



# Surgical Implant & Sinus Lift Kit

Features a compilation of clinicians' most popular procedural set-ups. Simply order a Collection today and receive high quality Hu-Friedy instruments and IMS cassettes together in procedural set-ups at significantly reduced prices. Create a professional, efficient working environment with the next generation of IMS Infinity cassettes.

This kit contains: IMS Infinity Cassette 1 DIN for 18

Iniskit Contains: IMS Infinity Cassette Fbin for Instruments, Ocean Blue, High Definition Mirror, Mirror Handle, Dressing Pliers No.18, Explorer, Probe. Tissue Pliers PermaSharb. Tissue Pliers

Handle

Surgical PermaSharp, Scalpel
No.5E, Buser Periosteal, Buser
Bone Scraper, Implant Probe
(7-10-13-15-18mm), Periosteal
Labanca, Lucas Surgical Curette
No.86, Molt curette No.2,
Jovanovic Plugger, Jovanovic
Implant Instrument, Surgical
Bone Curette Labanca, Palti
Sinus Lift Instrument Small,
Palti Sinus Lift Instrument,
Scissors Super Cut, NeedleHolder Simion, Cheek

DDH000

Retractor Simion.



# Swiss Perio Surgical Kit

This specialised surgical kit provides precision micro-instruments specifically designed for today's minimally invasive soft tissue grafting procedures. Surgical applications include soft tissue augmentation, root coverage, vestibular extension, stabilisation of mobile marginal tissue and ridge augmentation. The Kit includes: 1 DIN Surgical cassette for 18 instruments, black, Swiss Perio Microsurgical, Needle Holder Steel, Plier Straight Steel, Pcissor Steel, Scissors Steel, 45 degree ,Swiss Perio Mod. P24G, Black Line, Swiss Perio Mod. P24GSIC, Black Line, Swiss Perio Mod. PPSCHLEE, Black Line, Swiss Perio Mod. PPSCHLEE, Black Line, Swiss Perio Mod. P24GSP, Black Line, 84 Lucas Double Ended Surgical Curette,

Black Line, YG7/8 Younger-Good, EverEdge, Smooth handle, Colorvue Probe Handle, Smooth Handle, 12 Colorvue Probe tips, 12 Replacements Tips, Mini Scalpel Blade / Mirror Handle Stainless Steel, 5 European Style Round Scalpel Handle, N.6 Simple Steam Mirror Handle, Smooth Handle, 17 Dressing Pliers, Q2N Colour-Coded Nabers Probe, Smooth Handle, European Simple Steam N.5 Double-Sided HD Mirror, Single.

**DDH015** 

### **Surgical Procedure**

Instruments



### **Periosteal**

### Prichard Periosteal Short #6 Satin Steel Handle

Fits in IMS cassette. For reflecting and retracting the mucoperiosteum after incisions of the gingival tissue. Satin Steel handles offer a larger diameter and help lighten the instrument grasp and reduce hand fatigue.

#### DEE198



For delicate tissue retraction. Spear-shaped end initiates the flap by reflecting the papilla in a sulcular incision. Small paddle-shaped end continues the flap reflection.

#### DEE651



# **Surgical Curettes**

# 84 Lucas 2.2mm. DSC457 Lucas 85 2 5mm DSC456 Lucas 86 BlackLine DSC459

# Periotomes

#### No. 2 Anterior, Black Line

PT2 has a straight blade for anterior use. Thin sharp blades used to facilitate the removal of teeth with minimal damage to the surrounding alveolar bone by severing the periodontal ligament. May be used wherever an extraction is indicated, but are particularly important when considering the placement of the dental implant.

#### **DEE677**



# Probe UNC 15 Probe

Colour-coded with accurate, smoothly finished black markings.

#### DPR665



### **Pliers**

20 Corn Suture Pliers

15cm / 6".
DPI230

### 1x2 Adson Perma Sharp Tissue Pliers

Perma-Sharp tissue forceps/pliers used to grasp and stabilize soft tissue during suturing. 12cm.



### **College 2 Dressing Pliers**

Made of heavy gauge heat-treated stainless steel to prevent tip bending and misalignment.



# **Needle Holders**





### **Castro Micro Surgical**

Straight Diamond Dusted. 18cm.





### **Cheek Retractors**

**University of Minnesota** 



Cawood-Minnesota

DRH001

### **Scalpel Handle**

European Round No.5



# Scalpel blade remover

Allows simple, safe removal of blade from scalpel handle. Usable for all blade shapes



SHH000

# **Surgical Procedure**

Instruments

Hu-Friedy



# **Blackline Range**

Performance engineered coating for a harder, smoother surface for optimal edge retention and enhanced lubricity. Distinct black finish for enriched contrast between the instrument and the surgical site and underlying tissue. Reduced light reflection afforded by a matte finished handle and black working ends. Unique, smooth, large diameter lightweight handle for maximum comfort, reduced hand fatigue and increased control. Handcrafted with Immunity Steel® for optimal strength and corrosion resistance.

### **Tunnelling Knife No.2**

For incisions that remove or recontour soft tissue; also used to excise interproximal tissue.

DKH022

### **Buser periosteal Black line**

For reflecting and retracting the mucoperiosteum. DEE673

### Pritchard periosteal PR-3 black line

For reflecting & retracting the mucoperiosteum. DEE671

#### Lucas

Surgical currette used for curettage, cyst removal and debridement of tooth sockets.

DSC459

### 9 Molt periosteal elevator, black line

For reflecting & retracting the mucoperiosteum. DEE672



### UNODENT









# Disposable Scalpels

Disposable for convenience. Sealed individually.

SSS010 - No. 10 - Box 10 SSS011 - No. 11 - Box 10 SSS015 - No. 15 - Box 10



# Red Sterile, Carbon Steel Scalpels

SBS101 - No. 10 - Box 100 SBS111 - No. 11 - Box 100 SBS116 - No. 11P - Box 100 SBS112 - No. 12 - Box 100 SBS105 - No. 15 - Box 100 SBS115 - No. 15 - Box 100 SBS117 - No. 15T - Box 100



# Red Sterile, Stainless Steel Scalpels

SBS264 - No. SM64 - Box 25 SBS265 - No. SM65 - Box 25 SBS267 - No. SM67 - Box 25 SBS268 - No. SM68 - Box 25 SBS269 - No. SM69 - Box 25



# Photographic Mirror Handle



# Pure Reflect Mouth Mirrors

Box 12.

XXO050 - Size 0 Ø 14mm XXO051 - Size 1 Ø 16mm XXO052 - Size 2 Ø 18mm XXO053 - Size 3 Ø 20mm XXO054 - Size 4 Ø 22mm XXO055 - Size 5 Ø 24mm

XXO056 - Size 6 Ø 26mm





### **Soft Tissue Punch**

The circular scalpel or mucosal punch for soft tissues is an instrument used to create a complete incision of the mucous membrane flap of the diameter of the positioned implant. Omnia offers scalpels with 3 different diameters; the three sizes allow it to be adapted to the main implants on sale. The disposable circular scalpel causes less trauma to the mucous membrane than the classic scalpel because the cut is limited to the size of the implant cap. Box 5.

SDO012 - 4.1mm SDO014 - 5.2mm



### curasan

# **Safescraper Micross**

Disposable instrument with stainless steel blade to harvest autologous bone

- Easy and fast extraction of cortical bone grafts
- · Minimal invasive "Tunnel Technique"
- Ready to use, sterile instrument

**DEN590** 



## **Surgical Procedure**

**Dressings** 



### **Dermabond**

### **Prineo Skin Closure System**

Provides a flexible microbial barrier with 99% protection in vitro for 72 hours against organisms commonly responsible for SSIs.

**CCS100** - 60cm

CCS103 - 22cm

### Skin Adhesive - High Viscosity

May be used in place of sutures or staples and lends increased strength through layers by forming its own flexible and sterile waterresistant protective coating.

CCS101 - ProPen APP6

CCS102 - 0 5ml Vials



#### Curaspoi

# Absorbable Haemostatic Gelatin Sponge

Gamma sterilised, single-use,  $10 \times 10 \times 10$  mm, individually blister packed. This form is suitable for many of the most frequent surgical procedures in particular: parenchymal and thoracic surgery.



NJC012 - Pack 60