Surgical Procedure





Surgical Instruments

Every Anaqa product is subjected to a rigorous four-step quality control process. This robust process is the reason why Anaqa are able to confidently offer industry-leading warranties on all their instruments.

- Slim beaks for better access and hold
- Stunning finish with high corrosion resistance
- Four stage QC process
- Optional titanium coating for premium finish and maximum corrosion resistance
- High quality medical grade steel
- Smooth joints for stress-free working
- Ergonomic, lightweight design
- Industry-leading warranties

Luxation Instrument 5ST

Black Titanium Nano Coating. DEE803

Periotome

•

Black Titanium Nano Coating Available in sizes 1, 2 & 3 DEE814-816

Millers Bone File No. 52

Black Titanium Nano Coating Straight cut bone file used for final smoothing of bone. DEE813

Periotome Kit

Used to sever the gingival fibres and periodontal ligament. Once a tooth is removed the Periotome can be used to mobilise the soft tissue. Minimises treatment of the alveolar bone and comes with 3 blades to facilitate access throughout the mouth. Straight blades for anterior teeth and premolars. Angled blade for posterior teeth.

Contains - One standard stainless steel handle for detachable tips, one periotome blade narrow, one periotome blade wide and one periotome blade angled all in standard stainless steel.

DEE818

Blade P1 Wide DEE820 Blade P2 Narrow DEE819 Blade P3 Angled DEE821 Handle DEE822





Micro Surgical Basic Kit

The Kit includes:

IMS Infinity Cassette 1 DIN for 18 instruments, Black, High Definition Mirror, Mirror Handle, UNC probe 1-15mm, Micro Surgical Scalpel Handle SinusLine, Micro Tissue Plier, Periosteal Micro Surgery, Periosteal Trombelli, Micro Surgical Scissors Cortellini SuperCut. Needle holder

Castroviejo Perma Sharp, Suture plier, Urban Cheek Retractor, Labanca Malleable Retractors, Set of 3.

DDH008

Surgical Implant & Sinus Lift Kit

Features a compilation of clinicians' most popular procedural set-ups. Simply order a Collection today and receive high quality Hu-Friedy instruments and IMS cassettes together in procedural set-ups at significantly reduced prices. Create a professional, efficient working environment with the next generation of IMS Infinity cassettes. This kit contains: IMS Infinity Cassette 1 DIN for 18 instruments, Ocean Blue, High Definition Mirror, Mirror Handle, Dressing Pliers No.18, Explorer, Probe, Tissue Pliers PermaSharp, Tissue Pliers Surgical PermaSharp, Scalpel Handle

COLOR DE LE COL

No.5E, Buser Periosteal, Buser Bone Scraper, Implant Probe (7-10-13-15-18mm), Periosteal Labanca, Lucas Surgical Curette No.86, Molt curette No.2, Jovanovic Plugger, Jovanovic Implant Instrument, Surgical Bone Curette Labanca, Palti Sinus Lift Instrument, Scissors Super Cut, Needle-Holder Simion, Cheek Retractor Simion.

> Popular Product

Swiss Perio Surgical Kit

This specialised surgical kit provides precision micro-instruments specifically designed for today's minimally invasive soft tissue grafting procedures. Surgical applications include soft tissue augmentation, root coverage, vestibular extension, stabilisation of mobile marginal tissue and ridge augmentation. The Kit includes: 1 DIN Surgical cassette for 18 instruments, black, Swiss Perio Microsurgical, Needle Holder Steel, Plier Straight Steel, Plier Straight Steel, Scissor Steel, Scissors Steel, 45 degree ,Swiss Perio Mod. P24G, Black Line, Swiss Perio Mod. P24GSIC, Black Line, Swiss Perio Mod. PH2, Black Line,Swiss Perio Mod. PPSCHLEE, Black Line ,Swiss Perio Mod. P24GSP, Black Line, 84 Lucas Double Ended Surgical Curette, Black Line, YG7/8 Younger-Good,

Black Line, YG7/8 Younger-Good, EverEdge, Smooth handle, Colorvue Probe Handle, Smooth Handle, 12 Colorvue Probe tips, 12 Replacements Tips, Mini Scalpel Blade / Mirror Handle Stainless Steel, 5 European Style Round Scalpel Handle, N.6 Simple Steam Mirror Handle, Smooth Handle, 17 Dressing Pliers, Q2N Colour-Coded Nabers Probe, Smooth Handle, European Simple Steam N.5 Double-Sided HD Mirror, Single.

DDH015



Surgical Procedure



Periosteal Prichard Periosteal Short #6 Satin Steel Handle

Fits in IMS cassette. For reflecting and retracting the mucoperiosteum after incisions of the gingival tissue. Satin Steel handles offer a larger diameter and help lighten the instrument grasp and reduce hand fatigue.

DEE198



Periosteal Satin Steel

For delicate tissue retraction. Spear-shaped end initiates the flap by reflecting the papilla in a sulcular incision. Small paddle-shaped end continues the flap reflection.

DEE651



Surgical Curettes

84 Lucas





Periotomes No. 2 Anterior, Black Line

PT2 has a straight blade for anterior use. Thin sharp blades used to facilitate the removal of teeth with minimal damage to the surrounding alveolar bone by severing the periodontal ligament. May be used wherever an extraction is indicated, but are particularly important when considering the placement of the dental implant.

DEE677



Probe UNC 15 Probe

Colour-coded with accurate, smoothly finished black markings.

DPR665



Pliers 20 Corn Suture Pliers 15cm / 6"

DPI230

1x2 Adson Perma Sharp Tissue Pliers

Perma-Sharp tissue forceps/pliers used to grasp and stabilize soft tissue during suturing. 12cm. DPI205



College 2 Dressing Pliers

Made of heavy gauge heat-treated stainless steel to prevent tip bending and misalignment.

DPI200

Ronguer Friedman

14cm / 5½."

DFO354

Needle Holders Castrovieio Straight. 18cm. DNH003 Crile-Wood 15cm. Popular Product **DNH004**

Castro Micro Surgical

Straight Diamond Dusted. 18cm. DNH020

Mallet 1 Mead, 19cm / 71/2". With nylon facings. DCI021

Cheek Retractors University of Minnesota





Scalpel Handle

European Round No.5

SHH020

Scalpel blade remover

Allows simple, safe removal of blade from scalpel handle. Usable for all blade shapes SHHOOO



Page 33

Surgical Procedure

Hu-Friedy



Blackline Range

Performance engineered coating for a harder, smoother surface for optimal edge retention and enhanced lubricity. Distinct black finish for enriched contrast between the instrument and the surgical site and underlying tissue. Reduced light reflection afforded by a matte finished handle and black working ends. Unique, smooth, large diameter lightweight handle for maximum comfort, reduced hand fatigue and increased control. Handcrafted with Immunity Steel® for optimal strength and corrosion resistance.

Tunnelling Knife No.2

For incisions that remove or recontour soft tissue; also used to excise interproximal tissue. DKH022

Buser periosteal Black line

For reflecting and retracting the mucoperiosteum. DEE673

Pritchard periosteal PR-3 black line For reflecting & retracting the mucoperiosteum. DEE671

Lucas

Surgical currette used for curettage, cyst removal and debridement of tooth sockets. DSC459

9 Molt periosteal elevator, black line For reflecting & retracting the mucoperiosteum. DEE672



UNODENT









Disposable Scalpels

Disposable for convenience. Sealed individually.

SSS010 - No. 10 - Box 10 SSS011 - No. 11 - Box 10 SSS015 - No. 15 - Box 10



Photographic Mirror Handle

with rotating holder. XXO069

Red Sterile, Carbon Steel Scalpels

SBS101 - No. 10 - Box 100 SBS111 - No. 11 - Box 100 SBS116 - No. 11P - Box 100 SBS112 - No. 12 - Box 100 SBS105 - No. 15 - Box 100 SBS115 - No. 15 - Box 100 SBS117 - No. 15T - Box 100

Red Sterile, Stainless Steel Scalpels

SBS264 - No. SM64 - Box 25 SBS265 - No. SM65 - Box 25 SBS267 - No. SM67 - Box 25 SBS268 - No. SM68 - Box 25 SBS269 - No. SM69 - Box 25

Pure Reflect Mouth Mirrors

CONTRAN .

Box 12.

M6

XX0050 - Size 0 Ø 14mm XX0051 - Size 1 Ø 16mm XX0052 - Size 2 Ø 18mm XX0053 - Size 3 Ø 20mm XX0054 - Size 4 Ø 22mm XX0055 - Size 5 Ø 24mm XX0056 - Size 6 Ø 26mm





Soft Tissue Punch

The circular scalpel or mucosal punch for soft tissues is an instrument used to create a complete incision of the mucous membrane flap of the diameter of the positioned implant. Omnia offers scalpels with 3 different diameters: the three sizes allow it to be adapted to the main implants on sale. The disposable circular scalpel causes less trauma to the mucous membrane than the classic scalpel because the cut is limited to the size of the implant cap. Box 5.

SDO012 - 4.1mm SDO014 - 5.2mm



curasan

Safescraper Micross

Disposable instrument with stainless steel blade to harvest autologous bone

- · Easy and fast extraction of cortical bone grafts
- Minimal invasive "Tunnel Technique"
- Ready to use,

sterile instrument

DEN590



Surgical Procedure

Dressings

FTHICON

a Johnson Johnson company

Dermabond

Prineo Skin Closure System

Provides a flexible microbial barrier with 99% protection in vitro for 72 hours against organisms commonly responsible for SSIs. CCS100 - 60cm

CCS103-22cm

Skin Adhesive - High Viscosity

May be used in place of sutures or staples and lends increased strength through lavers by forming its own flexible and sterile waterresistant protective coating.

CCS101 - ProPen APP6

ccs102 - 0 5ml Vials



Absorbable Haemostatic **Gelatin Sponge**

Gamma sterilised, single-use, 10 x 10 x 10mm, individually blister packed. This form is suitable for many of the most frequent surgical procedures in particular: parenchymal and thoracic surgery.





